

Product Datasheet

ORF1 antibody (orb824362)

Description	ORF1 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rabbit hemorrhagic disease virus (strain Rabbit/Germany/FRG/1989) (Ra/LV/RHDV/GH/1989/GE) (RHDV-FRG) ORF1 protein
Target	ORF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27410
Expiration Date	12 months from date of receipt.