

Product Datasheet

PPNP antibody (orb824483)

Description	PPNP antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) ppnP protein
Target	ppnP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7LMJ4
Expiration Date	12 months from date of receipt.