

Product Datasheet

PYRC antibody (orb825070)

Description	PYRC antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks) pyrC protein
Target	pyrC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B1IV40
Expiration Date	12 months from date of receipt.