


Product Datasheet

Human CX3CL1 protein (orb82612)

Description	Fractalkine, also named neurotactin, is a novel chemokine recently identified through bioinformatics...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	CX3CL1
Preservatives	Lyophilized from a 0.2 µm filtered concentrated (1.0mg/ml) solution in 20mM PB, pH 7.4, 50mM NaCl. We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml. Further dilutions should be made in appropriate buffered solutions.
Protein Sequence	QHHGVTKCNI TCSKMTSKIP VALLIHYQQN QASCGKRAII LETQRHRLFC ADPKEQWVKD AMQHLDQRQA ALTRNG
Purity	> 97% pure by SDS-PAGE and HPLC analyses.
MW	8.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human CX3CL1 protei...