

Product Datasheet

YI9A antibody (orb830393)

Description	YI9A antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli Insertion element IS1397 uncharacterized 20.1 kDa protein
Target	Insertion element IS1397 uncharacterized 20.1 kDa
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q47309
Expiration Date	12 months from date of receipt.