

Product Datasheet

YAC1 antibody (orb830532)

Description	YAC1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Methanosaeta concilii (Methanotherx soehngeni) Uncharacterized 16.9 kDa protein in acs 3'region protein
Target	Uncharacterized 16.9 kDa protein in acs 3'region
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27096
Expiration Date	12 months from date of receipt.