

Product Datasheet

TCTP antibody (orb833411)

Description	TCTP antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Elaeis guineensis</i> var. <i>tenera</i> (Oil palm) TCTP protein
Target	TCTP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5J907
Expiration Date	12 months from date of receipt.