

Product Datasheet

ORF3b antibody (orb837042)

Description	ORF3b antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Tobacco streak virus (strain WC) (TSV) ORF3b protein
Target	ORF3b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03598
Expiration Date	12 months from date of receipt.