

Product Datasheet

ndhF antibody (orb837101)

| | |
|----------------------------|--|
| Description | ndhF antibody |
| Species/Host | Rabbit |
| Reactivity | Plant |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Gossypium barbadense (Sea-island cotton) (Egyptian cotton) ndhF protein |
| Target | ndhF |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | A0ZZ82 |
| Expiration Date | 12 months from date of receipt. |