

Product Datasheet

rpl23-A; antibody (orb837184)

Description	rpl23-A; antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Gossypium barbadense (Sea-island cotton) (Egyptian cotton) rpl23-A protein
Target	rpl23-A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0ZZ77
Expiration Date	12 months from date of receipt.