

Product Datasheet

G6PD antibody (orb837479)

Description	G6PD antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zea mays (Maize) Glucose-6-phosphate 1-dehydrogenase, cytoplasmic isoform protein
Target	Glucose-6-phosphate 1-dehydrogenase, cytoplasmic isoform
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P80619
Expiration Date	12 months from date of receipt.