

Product Datasheet

YAC9 antibody (orb837847)

Description	YAC9 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zea mays (Maize) Transposable element activator uncharacterized 23 kDa protein
Target	Transposable element activator uncharacterized 23 kDa
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03936
Expiration Date	12 months from date of receipt.