

Product Datasheet

UC31 antibody (orb837974)

Description	UC31 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zea mays (Maize) Unknown protein from spot 67 of 2D-PAGE of etiolated coleoptile protein
Target	Unknown protein from spot 67 of 2D-PAGE of etiolated coleoptile
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P80637
Expiration Date	12 months from date of receipt.