

Product Datasheet

YMS4 antibody (orb838004)

Description	YMS4 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zea mays (Maize) Uncharacterized 29 kDa protein in mitochondrial S-1 DNA protein
Target	Uncharacterized 29 kDa protein in mitochondrial S-1 DNA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10580
Expiration Date	12 months from date of receipt.