

## Product Datasheet

### Os10g0138100 antibody (orb838247)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Os10g0138100 antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Plant   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, WB   |
| <b>Immunogen</b>           | Recombinant <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) Os10g0138100 protein  |
| <b>Target</b>              | Os10g0138100  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol  |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q7XH05</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |