

Product Datasheet

GAPC1 antibody (orb839748)

Description	GAPC1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) GAPC1 protein
Target	GAPC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q0J8A4
Expiration Date	12 months from date of receipt.