

Product Datasheet

HSFA2D antibody (orb840100)

Description	HSFA2D antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Oryza sativa</i> subsp. japonica (Rice) HSFA2D protein
Target	HSFA2D
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8H7Y6
Expiration Date	12 months from date of receipt.