

Product Datasheet

P9 antibody (orb840427)

Description	P9 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rice dwarf virus (isolate Fujian) (RDV) Minor outer capsid protein P9 protein
Target	Minor outer capsid protein P9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q85440
Expiration Date	12 months from date of receipt.