

Product Datasheet

Os10g0467600 antibody (orb840602)

Description	Os10g0467600 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) Os10g0467600 protein
Target	Os10g0467600
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48600
Expiration Date	12 months from date of receipt.