

Product Datasheet

Os06g0617800 antibody (orb840726)

Description	Os06g0617800 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Oryza sativa</i> subsp. japonica (Rice) Os06g0617800 protein
Target	Os06g0617800
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q69XQ6
Expiration Date	12 months from date of receipt.