

Product Datasheet

Os02g0639900 antibody (orb842704)

Description	Os02g0639900 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Oryza sativa</i> subsp. japonica (Rice) Os02g0639900 protein
Target	Os02g0639900
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6H7E4
Expiration Date	12 months from date of receipt.