

Product Datasheet

Protein 6 6 antibody (orb842791)

Description	Protein 6 6 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rice yellow stunt virus (RYSV) (Rice transitory yellowing virus) 6 protein
Target	6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O70791
Expiration Date	12 months from date of receipt.