

Product Datasheet

p3 antibody (orb842793)

Description	p3 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rice hoja blanca virus (strain cr) (RHBV) (Rice hoja blanca virus (strain Costa Rica)) p3 protein
Target	p3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q67897
Expiration Date	12 months from date of receipt.