

Product Datasheet

gag-pol antibody (orb842949)

Description	gag-pol antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype D (isolate Z2/CDC-Z34) (HIV-1) gag-pol protein
Target	gag-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12499
Expiration Date	12 months from date of receipt.