

Product Datasheet

FOS-FOX antibody (orb843198)

Description	FOS-FOX antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant FBR murine osteosarcoma virus (FBR-MSV) (Finkel-Biskis-Reilly murine osteosarcoma virus) FOS-FOX protein
Target	FOS-FOX
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29176
Expiration Date	12 months from date of receipt.