

Product Datasheet

gag-pro-pol antibody (orb843431)

Description	gag-pro-pol antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human T-cell leukemia virus 1 (isolate Melanesia mel5 subtype C) (HTLV-1) gag-pro-pol protein
Target	gag-pro-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C211
Expiration Date	12 months from date of receipt.