

Product Datasheet

S1 antibody (orb843648)

Description	S1 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human spumaretrovirus (SFVcpz(hu)) (Human foamy virus) Putative uncharacterized protein S1 protein
Target	Putative uncharacterized protein S1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14352
Expiration Date	12 months from date of receipt.