

Product Datasheet

YPX2 antibody (orb843650)

Description	YPX2 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bovine leukemia virus (isolate Japanese BLV-1) (BLV) Putative uncharacterized protein PXBL-II protein
Target	Putative uncharacterized protein PXBL-II
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03413
Expiration Date	12 months from date of receipt.