

Product Datasheet

VPU antibody (orb843850)

Description	VPU antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype F1 (isolate 93BR020) (HIV-1) vpu protein
Target	vpu
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O70891
Expiration Date	12 months from date of receipt.