

Product Datasheet

GIP1 antibody (orb846356)

| | |
|----------------------------|--|
| Description | GIP1 antibody |
| Species/Host | Rabbit |
| Reactivity | Bacteria |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) GIP1 protein |
| Target | GIP1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P38229 |
| Expiration Date | 12 months from date of receipt. |