

Product Datasheet

GIP1 antibody (orb846367)

Description	GIP1 antibody
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Gibberella zeae</i> (strain PH-1 / ATCC MYA-4620 / FGSC 9075 / NRRL 31084) (Wheat head blight fungus) (<i>Fusarium graminearum</i>) GIP1 protein
Target	GIP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	I1RF62
Expiration Date	12 months from date of receipt.