

Product Datasheet

RBT4 antibody (orb849058)

Description	RBT4 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) RBT4 protein
Target	RBT4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5AB48
Expiration Date	12 months from date of receipt.