

Product Datasheet

RIO1 antibody (orb849128)

Description	RIO1 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite) RIO1 protein
Target	RIO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8SVI7
Expiration Date	12 months from date of receipt.