

---

## Product Datasheet

**4EBP1 (phospho Thr37/46) antibody (FITC)  
(orb8493)**

**Description**

Rabbit polyclonal to 4EBP1 (phospho

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human 4EBP1 around the phosphorylation site of Thr37/46
<b>Target</b>	EIF4EBP1
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	13kDa
<b>Uniprot ID</b>	<a href="#">Q13541</a>
<b>Dilution Range</b>	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
<b>Expiration Date</b>	12 months from date of receipt.