


Product Datasheet

Rheb protein (orb85053)

Description	Recombinant human Rheb expressed in E.coli as fusion protein with GST. The eukaryotic translational ...
Form/Appearance	liquid (Supplied in 25 mM HEPES pH 7.4, 5 mM MgCl 2and 10 μM GTP)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 47 kDa
Expiration Date	6 months from date of receipt.



Western blot analysis of using Rheb prot...