

## Product Datasheet

### HIV-1 p17/p24 protein (orb85296)

<b>Description</b>	The protein contains the full-length sequence of HIV-1 core proteins: Core protein p17 (matrix prote...
<b>Form/Appearance</b>	liquid (Supplied in PBS pH 7.8, NaCl, DTT and urea)
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	May be used in ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 95% (SDS-PAGE, RP-HPLC)
<b>Expiration Date</b>	6 months from date of receipt.