

Product Datasheet

HIV-1 p17/ p24/ gp120 protein (orb85297)

Description	The protein contains the full-length sequence of HIV-1 core and envelope proteins: Core protein p17 ...
Form/Appearance	liquid (Supplied in 25 mM Tris-HCl pH 8.0, 5 mM β -mercaptoethanol and 8 M urea)
Storage	store at -20 °C
Note	For research use only
Application notes	May be used in ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, RP-HPLC)
Expiration Date	6 months from date of receipt.