

## Product Datasheet

### TY1B-MR2 antibody (orb852974)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | TY1B-MR2 antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Bacteria  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, WB   |
| <b>Immunogen</b>           | Recombinant <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) TY1B-MR2 protein                               |
| <b>Target</b>              | TY1B-MR2  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q04670</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |