

Product Datasheet

HIV-1 p17/ p24/ gp120/ gp41 protein (orb85298)

Description	The protein contains the full-length sequence of HIV-1 core and envelope proteins: Core protein p17 ...
Form/Appearance	liquid (Supplied in 20 mM PBS pH 7.8, 20 mM NaCl, 1 mM DTT and 8 M urea)
Storage	store at -20 °C
Note	For research use only
Application notes	May be used in ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, HPLC)
Expiration Date	6 months from date of receipt.