

Product Datasheet

HIV-1 p24 core protein (orb85303)

Description	HIV-1 p24 Core is a non-glycosylated polypeptide chain, containing the full length sequence of the H...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl pH 8.5 and 8 M urea)
Storage	store at -20 °C
Note	For research use only
Application notes	May be used in ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, HPLC)
MW	Theoretical MW: 24 kDa
Expiration Date	6 months from date of receipt.