

## Product Datasheet

### EGT1 antibody (orb854065)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | EGT1 antibody                                                                                                                              |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Bacteria                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                  |
| <b>Immunogen</b>           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) egt1 protein                                               |
| <b>Target</b>              | egt1                                                                                                                                       |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                                                                                |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                     |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">O94632</a>                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |