

Product Datasheet

SPAC4H3.07c antibody (orb855756)

| | |
|----------------------------|---|
| Description | SPAC4H3.07c antibody |
| Species/Host | Rabbit |
| Reactivity | Bacteria |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC4H3.07c protein |
| Target | SPAC4H3.07c |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q10215 |
| Expiration Date | 12 months from date of receipt. |