

Product Datasheet

RMT2 antibody (orb855890)

Description	RMT2 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) rmt2 protein
Target	rmt2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q10170
Expiration Date	12 months from date of receipt.