

Product Datasheet

G6PDC antibody (orb858355)

Description	G6PDC antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Solanum tuberosum (Potato) Glucose-6-phosphate 1-dehydrogenase, chloroplastic protein
Target	Glucose-6-phosphate 1-dehydrogenase, chloroplastic
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q43839
Expiration Date	12 months from date of receipt.