

Product Datasheet

ATP7 antibody (orb858420)

Description	ATP7 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Solanum tuberosum (Potato) ATP synthase 27 kDa subunit, mitochondrial protein
Target	ATP synthase 27 kDa subunit, mitochondrial
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P80496
Expiration Date	12 months from date of receipt.