

Product Datasheet

PRS6B antibody (orb858635)

Description	PRS6B antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Solanum tuberosum (Potato) 26S proteasome regulatory subunit 6B homolog protein
Target	26S proteasome regulatory subunit 6B homolog
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54778
Expiration Date	12 months from date of receipt.