

Product Datasheet

hdgfl2 antibody (orb859822)

Description	Rabbit polyclonal antibody to hdgfl2. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) hdgfl2 protein
Target	hdgfl2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5XXA7
Expiration Date	12 months from date of receipt.