


## Product Datasheet

### Insulin protein (Active) (orb85999)

<b>Description</b>	Insulin protein (Active)
<b>Species/Host</b>	Yeast
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Target</b>	Insulin
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
<b>Storage</b>	Store at -20°C or - 80°C storage. Avoid repeated freeze / thaw cycles
<b>Note</b>	For research use only
<b>Purity</b>	98% , Affinity chromatography
<b>MW</b>	14.3 kDa
<b>Expiration Date</b>	6 months from date of receipt.



Western blot analysis of Insulin protein