

Product Datasheet

LEP-B antibody (orb860354)

Description	Rabbit polyclonal antibody to LEP-B.All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Oryzias latipes (Japanese rice fish) (Japanese killifish) LEP-B protein
Target	LEP-B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B9X0L5
Expiration Date	12 months from date of receipt.