

Product Datasheet

LRFN1 antibody (orb860367)

Description	Rabbit polyclonal antibody to LRFN1. All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) lrfn1 protein
Target	lrfn1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PGX3
Expiration Date	12 months from date of receipt.